

Manufactured by Careva

Crossit

A universal posture belt system.



Crossit will help achieve an upright sitting position which is required in order for the regular safety belt in vehicles to work properly.

Crossit can be used in any seat, provided there is an opening between the seat and the backrest. It is well suited for both private users and for public transportation.

If the passenger cannot or doesn't wish to have the belts crossed over the chest, you can strap the belts vertically and use the chest buckle to keep them fixed over the body.

Quick facts

- Universal - one size fits all
- Portable
- Requires no changes to the seat
- Easy to install
- Fits all seats with a space between the seat base and the backrest
- Can be used with wheelchairs with a stable backrest

To guarantee your safety Autoadapt products are designed and tested according to current directives and standards.

- 102581** Crossit
- 102582** Seat pad (option)
- 102583** Chest buckle (option)

NOTE! Crossit does NOT replace the vehicle safety belt.



Manufactured by Careva

Careva

A highly configurable posture belt system.



Careva will help achieve an upright sitting position which is required in order for the regular safety belt in vehicles to work properly. The system consists of range of single belts that can be combined as required for users of all ages who need additional support for an upright sitting position.

Careva can be used in any seat, provided there is an opening between the seat and the backrest.

Quick facts

- Highly configurable
- Available in Small, Medium and Large
- Requires no changes to the seat
- Easy to install
- Fits all seats with a space between the seat base and the backrest

To guarantee your safety Autoadapt products are designed and tested according to current directives and standards.

NOTE! Careva does NOT replace the vehicle safety belt.





Careva

Dimensions and models



The spinal belt is the foundation for all the Careva belt solutions and should always be installed first, before adding the additional belts.

- 101373** Careva Spinal belt, 3 channels, 2 buckles, Small
- 100955** Careva Spinal belt, 5 channels, 2 buckles, Medium/Large



The Combi kit is used to support the upper body and prevents the lower body from sliding forwards. The kit includes a spinal belt.

- 101670** Careva Combi kit, small
- 101671** Careva Combi kit, medium
- 101672** Careva Combi kit, large



The chest belt is used to position the upper body and prevents the person falling forwards or to the side.

- 100959** Careva chest (thoracic) belt, 80cm
- 100958** Careva chest (thoracic) belt, 100cm
- 100957** Careva chest (thoracic) belt, 120cm
- 100590** Careva chest (thoracic) belt, 135cm



The pelvic belt is used to prevent people from sliding forwards when sitting, or to keep very active passengers in place during travel.

- 100968** Careva pelvic T-belt, small
- 100587** Careva pelvic T-belt, medium
- 100969** Careva pelvic T-belt, large